ሐ	HAWAII HEALTH SYSTEMS	Department: Information Technology Division - Telemedicine	Policy No.: TEL 0018
€	C O R P O R A T I O N "Touching Lives Everyday"		Revision No.: N/A
	Policies and Procedures	Issued by: Dennis Sato Vice President & CIO	Effective Date: April 24, 2001
Subject: Changing Your Password in		Approved by:	Supersedes Policy: N/A
the Medreach Telemedicine Software		Thomas M. Driskill, Jr. President & CEO	Page: 1 of 2

I. **PURPOSE:** To provide instruction on how to change the user password in the Medreach Software program.

II. DEFINITIONS:

- **A.** User A person who has authority to utilize the Medreach software for Telemedicine purposes.
- **B.** Password A secret combination of alpha/numeric that is used to authenticate a user of the Telemedicine system.

III. PROCEDURE:

- **A.** Every user of the Medreach software will be issued a temporary password when they first come for training. That temporary password is: 'test'.
- **B.** After the user logs into the Medreach software for the first time, the temporary password must be changed.

C. Any time the user feels that his/her password has been compromised, the password should be changed as well.

- **D.** The process for changing the password is as follows:
 - 1. Log into Medreach using your user name and password.
 - 2. After getting into the program, scroll down on the left side of the screen until you reach an icon that says "Configuration."
 - 3. Single click on the icon.
 - 4. Single click on "User Information."
 - 5. You will see your user information.
 - 6. Single click on the purple 'edit' button at the bottom of the screen.
 - 7. Single click on the purple 'change password' button at the top portion of the screen.
 - 8. Enter the old password 'test.'

- 9. Enter your new password in both the 'Enter New Password' box and the 'Confirm New Password' box.
- 10. Single click the 'OK' button.
- 11. Single click the 'Save' button at the bottom of the screen.
- 12. Single click the 'Home' icon to return to the Telemedicine software.